

# AUTOMATIC BED HBD1K-A



# AUTOMATIC BED HBD1K-A

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD1K-A AUTOMATIC BED

4 pcs of European style guardrail, upside down structure

Luxury central control castors, brake the bed with one stamp

With 4 bumpers



## SPECIFICATIONS

Model	HBD1K-A
Back-rest lifting angle	0°-65°±5°
Foot-rest lifting angle	0°-35°±5°
Trendelenburg	12°
Reverse trendelenburg	12°
Dimension	2200x1020x420-780 mm

# Centrifugen

**Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)