

SEROLOGICAL PIPETTE PIC71-0300B



SEROLOGICAL PIPETTE PIC71-0300B

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

PIC71-0300B SEROLOGICAL PIPETTE

Easy to use

Medical grade original Polystyrene


RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation



SPECIFICATIONS

Model	PIC71-0300B
Volume	100 ml
Precision	2 %
Graduation	1
Color	Black
Sterile	Gamma
Package	Bulk pack
Pcs/bag	1 pc/bag
Pcs/ctn	300 pcs/carton
Unit	1

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com