

# SEROLOGICAL PIPETTE PIC71-1200C



# SEROLOGICAL PIPETTE PIC71-1200C

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

## PIC71-1200C SEROLOGICAL PIPETTE

Easy to use

Medical grade original Polystyrene

RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation



## SPECIFICATIONS

| Model                    | PIC71-1200C     |
|--------------------------|-----------------|
| Volume                   | 2 ml            |
| Precision                | 2 %             |
| Negative Graduation (ml) | -0.4            |
| Graduation               | 1 / 100         |
| Color                    | Green           |
| Sterile                  | Gamma           |
| Package                  | Bulk pack       |
| Pcs/bag                  | 25 pcs/bag      |
| Pcs/ctn                  | 1200 pcs/carton |
| Unit                     | 1               |

# Centrifugen

**Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)