

AUTO COAGULATION ANALYZER ANA13-140



AUTO COAGULATION ANALYZER ANA13-140

Auto Coagulation Analyzer used to test clotting items, Chromogenic, Immunologic measuring methods
Used in Clinical diagnostic, Laboratory, Pharmaceutical.
Also known as coagulation analyzer.

ANA13-140 AUTO COAGULATION ANALYZER



SPECIFICATIONS

| | |
|------------------------|--------------------|
| Model | ANA13-140 |
| Throughput | 140T/H |
| Analysis Method | Coagulation Method |
| Analysis Items | PT, APTT, TT, FIB |
| Sample Positions | 5 |
| Tilt reagent Positions | 6 |

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com