

# HORIZONTAL AUTOCLAVE AUT22-0146SA



# HORIZONTAL AUTOCLAVE AUT22-0146SA

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## AUT22-0146SA HORIZONTAL AUTOCLAVE



Product Image Coming Soon

### SPECIFICATIONS

Model	AUT22-0146SA
Capacity	146 L
Chamber Size (Wx H x L)	420x420x830 mm
External Size (WxHxL)	1290x1680x1100 mm
Net Weight	510 Kg
Power Consumption	2.5 KW
Power Supply	380 V 50 Hz

# Centrifugen

**Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)