

HORIZONTAL AUTOCLAVE AUT22-1500DB




HORIZONTAL AUTOCLAVE AUT22-1500DB

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

AUT22-1500DB HORIZONTAL AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	AUT22-1500DB
Capacity	1500 L
Chamber Size (Wx H x L)	680x1180x1870 mm
External Size (WxHxL)	1450 x1990x2080 mm
Net Weight	2111 Kg
Power Consumption	5 KW
Power Supply	380 V 50 Hz

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com