

MICROCENTRIFUGE CEN43-060



MICROCENTRIFUGE CEN43-060

This compact, powerful micro centrifuge is useful to separate or sediment small volume samples at high speed. The short-spin key, soft brake, FastTemp features optimizes the overall function of the centrifuges. Faults are diagnosed automatically and displays fault code. These hybrid micro-centrifuge can run PCR tubes, Capillary tubes along with standard tubes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover. Also known as Benchtop Centrifuge, Laboratory Centrifuge.

CEN43-060 MICROCENTRIFUGE

Quick release of rotor for easy rotor exchange

Gentle braking system for quick sample processing



SPECIFICATIONS

| | |
|---------------------------------------|-------------------|
| Model | CEN43-060 |
| Maximum Capacity (No of tubes x Vol.) | 8x1.5 ml / 2.0 ml |
| Maximum Speed RPM | 4000 & 6000 |
| Maximum RCF | 1200xg & 2000xg |
| Overall Dimension | 135x147x120 mm |
| Weight (Net/Gross) | 0.7 kg |
| Power Supply | AC100-240V |

OPTIONAL ACCESSORIES

| Accessory Code | Name | Description | Capacity |
|----------------|-------------------|-------------|----------------|
| LS55142 | Fixed-Angle Rotor | | 8x1.5/2.0ml |
| LS55143 | Strip Rotor | | 4-8x0.2ml |
| LS55144 | Adapter | | 0.5-0.6ml tube |
| LS55145 | Adapter | | 0.2ml tube |

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com