

AUTOMATIC BED HBD1L-A1



AUTOMATIC BED HBD1L-A1

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD1L-A1 AUTOMATIC BED

With the nurse controller and CPR function

Bed board adopts HPL board, x-ray photograph for whole body is available

2 step adjustment at shank section



SPECIFICATIONS

| | |
|-------------------------|----------------------|
| Model | HBD1L-A1 |
| Back-rest lifting angle | 0°-65°±5° |
| Foot-rest lifting angle | 0°-35°±5° |
| Bed panel extention | Manual/180mm |
| Load Bearing | 240 kg |
| Patient surface | 1950x910 mm |
| Dimension | 2170x1080x400~750 mm |

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com