

MANUAL TRANSFER STRETCHER HBD8F- MS



MANUAL TRANSFER STRETCHER HBD8F-MS

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

HBD8F-MS MANUAL TRANSFER STRETCHER

Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down Bright arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easy control the stretcher direction

Standard accessories: Mattress, I.V.Pole



SPECIFICATIONS

Model	HBD8F-MS
Back-rest lifting angle	0-60°
Height adjustment	860-560 mm
IV pole	Yes
Net Weight	82 kg

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com