

CUBE ICE MAKER

ICM21-0500



CUBE ICE MAKER ICM21-0500

Ice-making speed is highly vigorous and it is energy efficient. Due to self contained storage bin, it becomes more spacious and uniform cube shaped ice are produced rapidly.

Used in Healthcare, Food industries, Refrigerated perishables, Biochemical and biological studies, Medical laboratories, Beverage industries, Pharmaceutical.

Also known as Ice Cuber Machine, Laboratory Cube Ice Maker.

ICM21-0500 CUBE ICE MAKER

High grade stainless steel SUS304 is used which is anticorrosive and durable. The compact structure is simple and space-saving.

Fully automatically computer controlled, convenient for use.

Automatic detection and protection in case of operation, water deficiency, ice full and malfunction.

Highly efficient fluoride-free compressor, energy saving and environment-friendly.

Suitable for Pubs, KTV and dancing halls, gyms, cold drinks shops, milk tea shops, cafes, convenience stores, restaurants and other entertainment venues.

Split-type structure, large ice bin, large ice storage capacity and strong applicability.



SPECIFICATIONS

Model	ICM21-0500
Capacity	500 kg/24 h
Storage Bin Capacity	280 kg
Refrigerant	R134a / R404a
Condenser Type	Water cooled
Machine Head	760x600x740 mm
Cabinet	760x800x980 mm
Packing Machine Head	845x705x950 mm
Packing Cabinet	800x850x1060 mm
Weight	123/133 kg
Power Supply	220 V, 60 Hz

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com