

GENERAL PURPOSE INCUBATOR INC81-080



GENERAL PURPOSE INCUBATOR INC81-080

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

INC81-080 GENERAL PURPOSE INCUBATOR

- LCD PID control system
- International brand cooling fan
- Stainless steel chamber
- Cavity preheating technology and duct cycling system



SPECIFICATIONS

Model	INC81-080
Capacity	80 L
Temperature Range	RT +5 ~ 65°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±0.5%
Internal Dimension (WxDxH)	480x480x350 mm
Overall Dimension (WxDxH)	680x630x600 mm
Shelves	2
Weight	42 kg
Power Supply	220 V 50 Hz

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com