

HOLE-TYPE OR LAMP OPL3B-12



HOLE-TYPE OR LAMP OPL3B-12

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

OPL3B-12 HOLE-TYPE OR LAMP



SPECIFICATIONS

Model	OPL3B-12
Illumination	50000 Lux
Depth of light beam	≥ 500 mm
The value of the total irradiance	$< 6\text{mW} / \text{m}^2\text{lx}$
Total irradiance	$\leq 600\text{W}/\text{m}^2$
Color rendering Index (Ra)	$100 \geq \text{Ra} \geq 85$
Color Temperature (K)	4000 ± 500 K
Power supply	$220 \pm 10\%$ 50 Hz DC24V
Power of bulb	125 W
Input power	150 VA

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com