

# MANUAL OR TABLE OPT11-2000



# MANUAL OR TABLE OPT11-2000

Manual Operation table are suitable for various operation tables conducted by medical unit. The OR table and accessories are ergonomic and can be used across multiple disciplines. Various actions are controlled by hydraulic pneumatic system. 360° rotation helps every surgical requirement.

Used in Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology, abdominal, gall bladder, plastic surgery and more.

Also known as Manual surgical table, Manual OR Table, Manual operating table.

## OPT11-2000 MANUAL OR TABLE

Various actions are controlled by hydraulic pneumatic system

All movements of the table are manipulated using the operating handle under the head plate

Imported Y type sealing ring is durable

The table surface adopts quality carbon plastic plate which is easily for conducting X-ray examination, equipped with coincidence rail

Both sides of the bed rail for the high quality 304 stainless steel

The base adopts imported carbon cloth which can be easily cleaned and sterilized

The mattress is wrapped with imported artificial leather. High-density memory sponge is formed in one piece, with no seam

Foot treadle brake device which is easily for moving operating table



## SPECIFICATIONS

Model	OPT11-2000
Table Length	2000±50 mm
Table Width	500±20 mm
Height Adjustment	(700~950) ±50 mm
Left tilt	≥ 15°
Right tilt	≥ 15°
Trendelenburg	≥ 20°
Reverse trendelenburg	≥ 15°
Head board upwards	≥ 30°
Head board downwards	≥ 90°
Backboard upwards	≥ 70°
Backboard downwards	≥ 15°
Waist board lifting	≥ 70°
Leg board downwards	≥ 90°
Leg board outwards	≥ 180°
Main Voltage AC	220V±22V, 50Hz±1Hz
Input Power VA	150±10%



***Centrifugen***

**Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)