

DRYING OVEN OVE11-040



DRYING OVEN OVE11-040

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

OVE11-040 DRYING OVEN

Precision LCD

Cooling fan

Stainless steel inner chamber

Intelligent PID control system

Best quality oven



SPECIFICATIONS

Model	OVE11-040
Capacity	40 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Temperature Uniformity	±1°C
Temperature Resolution	±0.5°C
Timer Range	0-999 min
Internal Dimension (WxDxH)	350x350x320 mm
Overall Dimension (WxDxH)	640x550x530 mm
Noise Level	≤55 db
Weight	36/39 kg
Power	1000 W
Power Supply	220 V, 50 Hz

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com