

DRYING OVEN OVE12-225



DRYING OVEN OVE12-225

Heating Oven, Laboratory Drying Oven

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory general Purpose Oven.

OVE12-225 DRYING OVEN

Micro-computer intelligent control system with LED display

Forced air convection maintains temperature uniformity

High quality stainless steel inner chamber with electrostatic spray treated exterior

Large observation window

Easy to clean semicircular corners with adjustable shelves

Silicon rubber sealing ring prevents heat loss

Over temperature alarm with timer



SPECIFICATIONS

Model	OVE12-225
Capacity	225 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±1°C
Temperature Uniformity	±7.0°C at 400°C
Temperature Resolution	±1°C
Timer Range	0-999 min
Internal Dimension (WxDxH)	750x600x500 mm
Overall Dimension (WxDxH)	980x965x770 mm
Package Size	1100x1045x850 mm
Shelves	2
Weight	93/102kg
Power	2500 W
Power Supply	220V, 60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS54389	Intelligent program temperature controller	
LS54401	Independent limit of temperature controller	

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com