

# THERMOSTATIC DRYING OVEN OVE13-070



# THERMOSTATIC DRYING OVEN OVE13-070

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE13-070 THERMOSTATIC DRYING OVEN

ALLSENSTM programmable PID control

ALLHEAT™ cavity warm-up technology

ALLFLOW™ perfect air current cycling

Pro-INSULATION™ Insulation technology

SUS304 stainless steel internal material

Build-in memory menu, real-time monitoring

Display parameters: Temperature, Air velocity, Run time, Run/Stop



## SPECIFICATIONS

Model	OVE13-070
Capacity	70 L
Temperature Range	RT+5 – 250°C
Ambient Temperature	10-30°C
Temperature Accuracy	0.1°C
Temperature Fluctuation	±1°C(50 - 240°C)
Temperature Uniformity	2%(50 - 240°C)
Ambient Humidity	<70%
Timer Range	0-99 h or 0-9999 min
Convection Mode	Forced Convection
Internal Dimension (WxDxH)	450x450x350 mm
Overall Dimension (WxDxH)	670x680x593 mm
Package Size	790x760x680 mm
Shelves	2
Weight	48/51 kg
Power	970 W
Power Supply	Single phase AC220V/50Hz



# ***Centrifugen***

## **Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)