

96 WELL GRADIENT THERMAL CYCLER PCR14-096C



96 WELL GRADIENT THERMAL CYCLER

PCR14-096C

Engineered by finest quality and leading edge technology according to the advance technology and market norms under the direction of competent experts. Simple, intuitive programming, cost-efficient, fast setup and convenient to use makes it an ideal choice.

Used in Analytical Laboratories, Molecular biology, Gene amplification, Gene Expression, Research, Development, Food Science, Pharmaceutical, Life Science, Animal Diagnostics.

Also known as Gradient PCR Machine, Laboratory Gradient PCR Thermal Cyclers, Gradient PCR Thermocycler, Laboratory Gradient PCR Thermocycler.

PCR14-096C 96 WELL GRADIENT THERMAL CYCLER

- Six pieces of long service life Peltier heating units and form 3 circuits to control 3 temperature zones
- Reinforced aluminum module with anodizing technology can keep rapid heating-conducting property and have enough corrosion resistance
- High heating and cooling rate, max. Ramping rate 4.5 °/s, can save your precious time
- Stepless adjustable hot lid, fit tubes of different heights to avoid tube melt and evaporation
- TFT color touch-screen with graphical display provides easy use for setting up and monitoring
- Built-in 11 standard program file template, can quickly edit the required files
- Folder management, user can build directory
- The running program and left time can be displayed in real time, allow to edit file when program is running
- One-click quick incubation function can meet experiment's needs such as denaturation, enzyme cutting/enzyme-link and ELISA
- Internal flash memory for 10000 typical PCR files in free configurable folders
- Hot lid temperature and hot lid work mode can be set to meet different experiment's need
- Automatic restart after power failure. When power is restored it can continue to run unfinished program
- GLP report records every step to provide accurate data support for experiment result analysis
- User Login Management, three-tier permission, password protection function to ensure data security
- Compatible with devices such as Mouse and Keyboard and capable to transfer data and perform software updates via USB Drive
- Support USB and LAN to update software
- One computer can control many sets of PCR via network connection
- Min type Bluetooth printer as an option, easily record information
- Support email-alert function when experiment is over



SPECIFICATIONS

Model	PCR14-096C
Sample Capacity	96×0.2 ml+77×0.5 ml

Temperature Range	0°C~99.9°C
Temperature Increment/Decrement	0.1~10.0°C
Hold at 4°C	Forever
Max. ramp rate	0.1°C~2.5°C
Max Heating Rate	4.5°C
Max Cooling Rate	4°C / s
Display Interface	LCD, 8", 800×600
Display Resolution	0.1°C
Uniformity	≤±0.3°C
Accuracy	≤±0.2°C
Thermal Gradient Accuracy	≤±0.3°C
Gradient Spread	1~30°C
Gradient Uniformity	≤±0.2°C
Hot Lid Temperature	30°C~110°C
Height of hot Lid	Stepless Adjustable
Max.No.of Cycle	100
Program Storage	10000+(USB Flash)
Max Program Steps	30
Communication	USB2.0 , LAN
Temp Control Mode	Block, tube
Time Increment/Decrement	1 sec ~600 sec
Pause Function	Yes
Auto Data Protection	Yes
Dimension (W×D×H)	270x390x255 mm
Power	600 W
Weight	9 kg
Power Supply	85~264 V AC , 47~63 Hz

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com