

STACKABLE 300ML TIP BOX PIC66-00300S



STACKABLE 300ML TIP BOX PIC66-00300S

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology. Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

PIC66-00300S STACKABLE 300ML TIP BOX

Made from medical grade Polypropylene
Slippy surface with excellent design
DNase-free
Rnase-free
Non-pyrogenic
Filter and non-filter for choice
Bag and rack are available to pack



SPECIFICATIONS

Model	PIC66-00300S
Volume	300 μ l
Unit	20

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com