

# 300ML TIP BOX WITH ULTRA LOW RETENTION TIP PIC6A-00300UL



# 300ML TIP BOX WITH ULTRA LOW RETENTION TIP PIC6A-00300UL

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology. Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

## PIC6A-00300UL 300ML TIP BOX WITH ULTRA LOW RETENTION TIP

- Made from medical grade Polypropylene
- Slippy surface with excellent design
- DNase-free
- Rnase-free
- Non-pyrogenic
- Filter and non-filter for choice
- Bag and rack are available to pack



### SPECIFICATIONS

Model	PIC6A-00300UL
Volume	300 µl
Unit	10



#### Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: info@centrifugen.com | Website: centrifugen.com