

SEROLOGICAL PIPETTE PIC71-0800C



SEROLOGICAL PIPETTE PIC71-0800C

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

PIC71-0800C SEROLOGICAL PIPETTE

Easy to use

Medical grade original Polystyrene

RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation



SPECIFICATIONS

| Model | PIC71-0800C |
|--------------------------|----------------|
| Volume | 5 ml |
| Precision | 2 % |
| Negative Graduation (ml) | -2.5 |
| Graduation | 1 / 10 |
| Color | Blue |
| Sterile | Gamma |
| Package | Individual |
| Pcs/bag | 25 pcs/bag |
| Pcs/ctn | 800 pcs/carton |
| Unit | 1 |

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com