

PORTABLE REFRIGERATOR REF32-340



PORTABLE REFRIGERATOR REF32-340

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

REF32-340 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



SPECIFICATIONS

Model	REF32-340
External Size (WxDxH)	340x230x280 mm
Internal Size (WxDxH)	226x140x210 mm
Capacity	6 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	24 h
Material	
External	PP
Internal	PU
Ice Vessels	6(CMB-400)
Max Loading	6 kg
Optional Accessory	LED display thermometer
Net Weight	1.9 kg
Gross Weight	5 kg
Package Size (WxDxH)	360x250x300 mm

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com