

# PORTABLE REFRIGERATOR REF32-615



# PORTABLE REFRIGERATOR REF32-615

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## REF32-615 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



## SPECIFICATIONS

Model	REF32-615
External Size (WxDxH)	615x395x450 mm
Internal Size (WxDxH)	525x300x355 mm
Capacity	54 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	24 h
Material	
External	PP
Internal	PU
Ice Vessels	4(CMB-1500)
Max Loading	11 kg
Optional Accessory	LED display thermometer
Net Weight	7.8 kg
Gross Weight	17 kg
Package Size (WxDxH)	800x550x510 mm

# Centrifugen

**Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)