

MAGNETIC STIRRER STR13-001



MAGNETIC STIRRER STR13-001

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR13-001 MAGNETIC STIRRER

Maintenance free brushless DC motor guarantees long service

Rare earth magnets give powerful stirring

Stepless speed regulation

Ultra silence and vibration-free design.

PET top provides excellent chemical resistance

Shaft can be fitted with a pH sensor or sensors for other purposes



SPECIFICATIONS

Model	STR13-001
Capacity	1.5 L
Stirring Positions	1
Speed Range	300-2000 prm
Permissible Ambient Temperature	5-40°C
Humidity	80%RH
Work Plate Size	Φ120 mm
Stirrer Dimension	40 mm
Platform Material	PET
Motor Type	Brushless DC motor
Protection Class	IP21
Overall Dimension	150Wx145Dx50H mm
Weight	0.3 kg
Power	7.2 W
Power Supply	100-120/200-240 V, 50/60 Hz

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@centrifugen.com | Website: centrifugen.com