

# BIOSAFETY TRANSPORT BOX ULT91-06



# BIOSAFETY TRANSPORT BOX ULT91-06

The Biosafety Transport Box is a non-electric cooler used for the transportation of biological products requiring physical insulation, such as biological samples, pathogenic microorganism (virus) species, blood, etc. Used in It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism (virus) species, blood, etc.

## ULT91-06 BIOSAFETY TRANSPORT BOX

Real time temperature value display

Refrigeration and temperature control to ensure the consistency of samples

Easy to operate and could be used repeatedly

Meet the basic requirements of cold chain transportation and storage of vaccines



## SPECIFICATIONS

Model	ULT91-06
Inner Capacity	6 L
Outer size	330x235x275 mm
Product Configuration	1 box 1 strap 1 digital thermometer 1 95K Pa bio-safety transport tank 1 sponge reagent tube holder 1 piece of absorbent cotton 4 crystal ice boxes
Material	PP+PE
Insulation Material	PU
Gross weight	2.8 kg
Packing size	380x260x290 mm

# Centrifugen

**Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)