

# ICE FLAKER



# ICE FLAKER

Suits a variety of applications including fisheries, laboratories and industry. Once formed, Ice flakes are dry, clean, beautiful in shape. Relatively big contact area and good mobility can ensure their full contact with refrigerated objects to realize good preservation effect. Great adaptability, sanitary and convenient.

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

## ICM31 ICE FLAKER

High quality stainless steel construction

Fully automatic operation system

Fluorine-free foam thermal cover ensures good insulation

High speed compact clear flake ice-production

Advanced self-checking mode function

Power-overloading, water shortage and ice storage full protection

Large output and energy efficient

Easy to clean storage bin



## SPECIFICATIONS

| Model                | ICM31-020      | ICM31-030  | ICM31-040      | ICM31-050  |
|----------------------|----------------|------------|----------------|------------|
| Capacity             | 20 kg/24 h     | 30 kg/24 h | 40 kg/24 h     | 50 kg/24 h |
| Storage Bin Capacity | 10 kg          |            | 15 kg          |            |
| Refrigerant          | R134a          |            |                |            |
| Condenser Type       | Air cooled     |            |                |            |
| Overall Dimension    | 330x470x605 mm |            | 400x543x720 mm |            |
| Packing Dimension    | 405x555x665 mm |            | 470x615x780 mm |            |
| Weight               | 28/31 kg       |            | 36/39 kg       |            |
| Power                | 280 W          | 300 W      | 360 W          | 380 W      |
| Power Supply         | 220 V, 60 Hz   |            |                |            |

## OPTIONAL ACCESSORIES

| Accessory Code | Name         | Description |
|----------------|--------------|-------------|
| LS54713        | Water filter |             |



ICM31-020



ICM31-030



ICM31-040



ICM31-050

## ICM32 ICE FLAKER

- High quality stainless steel construction
- Fully automatic operation system
- Fluorine-free foam thermal cover ensures good insulation
- High speed compact clear flake ice-production
- Advanced self-checking mode function
- Power-overloading, water shortage and ice storage full protection
- Large output and energy efficient
- Easy to clean storage bin



### SPECIFICATIONS

| Model                | ICM32-060      | ICM32-070  | ICM32-085  |
|----------------------|----------------|------------|------------|
| Capacity             | 60 kg/24 h     | 70 kg/24 h | 85 kg/24 h |
| Storage Bin Capacity | 25 kg          |            |            |
| Refrigerant          | R134a          |            |            |
| Condenser Type       | Air cooled     |            |            |
| Overall Dimension    | 400x510x825 mm |            |            |
| Packing Dimension    | 470x590x875 mm |            |            |
| Weight               | 42/46 kg       |            | 45/49 kg   |
| Power                | 420 W          | 460 W      | 480 W      |
| Power Supply         | 220 V, 60 Hz   |            |            |

| Model                | ICM32-100       | ICM32-120   | ICM32-200       |
|----------------------|-----------------|-------------|-----------------|
| Capacity             | 100 kg/24 h     | 120 kg/24 h | 200 kg/24 h     |
| Storage Bin Capacity | 40 kg           |             | 60 kg           |
| Refrigerant          | R134a           |             |                 |
| Condenser Type       | Air cooled      |             |                 |
| Overall Dimension    | 500x611x940 mm  |             | 680x611x1100 mm |
| Packing Dimension    | 570x685x1000 mm |             | 750x675x1150 mm |
| Weight               | 58/63 kg        |             | 80/85 kg        |
| Power                | 520 W           | 550 W       | 1040 W          |
| Power Supply         | 220 V, 60 Hz    |             |                 |

## OPTIONAL ACCESSORIES

| Accessory Code | Name         | Description | For Models                                 |
|----------------|--------------|-------------|--|
| LS54713        | Water filter |             | ICM32-060, ICM32-070, ICM32-085, ICM32-100 |



ICM32-060



ICM32-070



ICM32-085



ICM32-100



ICM32-120



ICM32-200

## ICM33-300 ICE FLAKER

- High quality stainless steel construction
- Fully automatic operation system
- Fluorine-free foam thermal cover ensures good insulation
- High speed compact clear flake ice-production
- Advanced self-checking mode function
- Power-overloading, water shortage and ice storage full protection
- Large output and energy efficient
- Easy to clean storage bin



## SPECIFICATIONS

|                      |                 |
|----------------------|-----------------|
| Model                | ICM33-300       |
| Capacity             | 300 kg/24 h     |
| Storage Bin Capacity | 150 kg          |
| Refrigerant          | R134a           |
| Condenser Type       | Air cooled      |
| Overall Dimension    | 560x800x1700 mm |
| Packing Dimension    | 650x850x1750 mm |
| Weight               | 95/105 kg       |

|              |              |
|--------------|--------------|
| Power        | 1360 W       |
| Power Supply | 220 V, 60 Hz |

# ***Centrifugen***

## **Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)