

WHEELCHAIR



WHEELCHAIR

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient Indoor-outdoor locomotion, Patient transfer.

Also known as Patient transfer wheelchair, Transport wheelchairs.

HBD91 MANUAL WHEELCHAIR

It is a type of steady and simple Manual Wheelchair



SPECIFICATIONS

Model	HBD91-MW0	HBD91-MW1	HBD91-MW2
Wheelchair length	970 mm		1020 mm
Opened wheelchair width	660 mm		630 mm
Wheelchair Height	920 mm		890 mm
Closed wheelchair width	350 mm		345 mm
Seat Width		430 mm	
Front wheel		6" inch	
Rear wheel		24" inch	
Seat Height	470 mm		460 mm
Seat depth		400 mm	
Backrest Height		450 mm	
Armrest Height		210 mm	
Distance between armrest	425 mm		445 mm
Height between footrests and ground	120 mm		80 mm
Net Weight	12.7 Kg	12.94 Kg	13.28 Kg

Model	HBD91-MW3	HBD91-MW4	HBD91-MW5
Wheelchair length	950 mm	855 mm	1250 mm
Opened wheelchair width	540 mm	535 mm	690 mm
Wheelchair Height	870 mm	855 mm	1210 mm
Closed wheelchair width	280 mm	250 mm	365 mm
Seat Width	405 mm		450 mm

Front wheel	6" inch		8" inch
Rear wheel	16" inch	14" inch	24" inch
Seat Height	445 mm	480 mm	520 mm
Seat depth	420 mm	430 mm	455 mm
Backrest Height	380 mm	350 mm	450 mm
Armrest Height	220 mm	190 mm	240 mm
Distance between armrest	435 mm	440 mm	450 mm
Height between footrests and ground	95 mm	60 mm	100 mm
Net Weight	10.24 Kg	13.75 Kg	25.4 Kg

Model	HBD91-MW6	HBD91-MW7	HBD91-MW8
Wheelchair length	1055 mm	830 mm	1180 mm
Opened wheelchair width	650 mm	580 mm	660 mm
Wheelchair Height	920 mm	780 mm	1320 mm
Closed wheelchair width	350 mm		360 mm
Seat Width	455 mm	360 mm	460 mm
Front wheel	8" inch	4" inch	8" inch
Rear wheel	24" inch		
Seat Height	510 mm	470 mm	520 mm
Seat depth	450 mm	380 mm	460 mm
Backrest Height	390 mm	330 mm	490 mm
Armrest Height	200 mm	215 mm	230 mm
Distance between armrest	450 mm	415 mm	450 mm
Height between footrests and ground	100 mm		115 mm
Net Weight	17.4 Kg	13 Kg	23.35 Kg



HBD91-MW0



HBD91-MW1



HBD91-MW2



HBD91-MW3



HBD91-MW4



HBD91-MW5



HBD91-MW6



HBD91-MW7



HBD91-MW8

HBD92 ELECTRIC WHEELCHAIR

It is a type of steady, simple and Practical Electric Wheelchair



SPECIFICATIONS

Model	HBD92-EW0	HBD92-EW1	HBD92-EW2
Wheelchair length	1055 mm	1050 mm	-
Opened wheelchair width	670 mm	685 mm	-
Wheelchair Height	880 mm	940 mm	-
Closed wheelchair width	500 mm	420 mm	-
Seat Width	445 mm		460 mm
Front wheel	8" inch		-
Rear wheel	22" inch	12" inch	-
Seat Height	480 mm	495 mm	460 mm
Seat depth	430 mm	420 mm	-
Backrest Height	400 mm	430 mm	-
Armrest Height	210 mm	205 mm	-
Distance between armrest	440 mm		-
Height between footrests and ground	100 mm		-
Km of Endurance	> 20 Km		
Gradeability (°)	> 6	-	-
Speed per Hour	> 6 Km/h		
Net Weight (Without battery)	29.1 Kg	34 Kg	-
Bearing Load	-	100 Kg	-
Gradeability	-	> 6°	-
Total Length	-	-	1000 mm
Total Width	-	-	620 mm
Total Height	-	-	890 mm

Seat Depth	-	-	400 mm
Maximum Load	-	-	160 Kg
Maximum Speed	-	-	7.2 Km/h
Mileage	-	-	17 Km
Damping System	-	-	2
Motor Power	-	-	24/250x2 V/W
Charger	-	-	24 (V2/A)
Battery	-	-	24/10 (V/Ah)
Charging Time	-	-	3-5 (h)
Front Wheel	-	-	8" inch
Rear Wheel	-	-	12" inch
Wheel Base	-	-	650 mm
Maximum Gradeability	-	-	12 °
Net Weight	-	-	25 Kg



HBD92-EW0



HBD92-EW1



HBD92-EW2

Centrifugen

Centrifugen

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@centrifugen.com | Website: centrifugen.com