

# STEREO ZOOM MICROSCOPE



# STEREO ZOOM MICROSCOPE

Stereo Zoom Microscope also known as the dissecting microscope is type of optical microscope that allows the user to see a three-dimensional view of a specimen.

Used in Pharmaceutical, University, Medical, Laboratory, Research.

## MSC41 STEREO ZOOM MICROSCOPE

Excellent optical function with infinite optical system.

Innovative stand structure and sharp image display.

Convenient operation and special design for viewing incubating cell tissue.

Valid mould proof design.

It can be used under the high temperature and high humidity environment.



Product Image Coming Soon

## SPECIFICATIONS

Model	MSC41-045M	MSC41-045MT	MSC41-020ST
Viewing Head			
Binocular view head inclined at 60°	Yes	No	-
Trinocular view head inclined at 45°	No	Yes	-
Digital view head inclined at 45°	No	Yes	-
Eyepiece			
High point wide field adjustable eyepiece WF10x22 mm	Yes		-
MSC41-020ST	-	-	-
MSC41-040ST	-	-	-
MSC41-060ST	-	-	-
Zoom Objective			
0.67 X—4.5 X	Yes		-
Zoom Ratio			
6.7:1	Yes		-
Auxiliary Objective			
0.5X /177 mm, 1.5X /47 mm, 2X/26 mm	Yes		-
Working Distance			
115 mm	Yes		-
MSC41-020ST	-	-	-
MSC41-040ST	-	-	-
MSC41-060ST	-	-	-
Focusing Bracket			
Focusing frame and adjustment knob	Yes		-
Interpupillary Distance			
55-75 mm / 52-75 mm	Yes		-

Up Illumination			
LED lamp	Yes	-	
Down Illumination			
LED lamp	Yes	-	
Power Supply			
AC110/220 V $\pm$ 10%, 50/60 Hz	Yes	-	
MSC41-020ST	-	-	-
MSC41-040ST	-	-	-
MSC41-060ST	-	-	-
Package			
Cartoon with foam	Yes	-	
Package Size	500x260x370 mm	300x200x150 mm	
Gross Weight	8 kg	3 kg	
Illumination	-	-	Natural light Binocular
View Head And Body	-	-	Straight Head
Objective	-	-	2X

Model	MSC41-040ST	MSC41-060ST
Viewing Head		
Binocular view head inclined at 60°	-	-
Trinocular view head inclined at 45°	-	-
Digital view head inclined at 45°	-	-
Eyepiece		
High point wide field adjustable eyepiece WF10x22 mm	-	-
MSC41-020ST	-	-
MSC41-040ST	-	-
MSC41-060ST	-	-
Zoom Objective		
0.67 X—4.5 X	-	-
Zoom Ratio		
6.7:1	-	-
Auxiliary Objective		
0.5X /177 mm, 1.5X /47 mm, 2X/26 mm	-	-
Working Distance		
115 mm	-	-
MSC41-020ST	-	-
MSC41-040ST	-	-
MSC41-060ST	-	-
Focusing Bracket		
Focusing frame and adjustment knob	-	-
Interpupillary Distance		
55-75 mm / 52-75 mm	-	-
Up Illumination		
LED lamp	-	-
Down Illumination		
LED lamp	-	-
Power Supply		
AC110/220 V $\pm$ 10%, 50/60 Hz	-	-

MSC41-020ST	-	-
MSC41-040ST	-	-
MSC41-060ST	-	-
Package		
Cartoon with foam	-	-
Package Size	360x170x290 mm	400x250x360 mm
Gross Weight	6 kg	8 kg
Illumination	Tranilluminator / inclined light	Natural light
View Head And Body	Inclined viewing Up-rotatable body	Inclined viewing Body 360°rotatable
Objective	2X&4X, 90°rotated	2X&4X, 360°rotated



Product Image Coming Soon

MSC41-045M



MSC41-045MT



MSC41-020ST



MSC41-040ST



MSC41-060ST

# Centrifugen

**Centrifugen**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@centrifugen.com](mailto:info@centrifugen.com) | Website: [centrifugen.com](http://centrifugen.com)